



MAKING SCIENCE
AFFORDABLE

RIGID FRAME 384-WELL SKIRTED PCR PLATE

Instrument Compatibility List

telephone: +44 (0)1279 940 983
fax: +44 (0)1707 240 451
email: support@3crbio.com
web: <https://3crbio.com/>

Unit 10, West Point Business Park, West Road, Harlow,
Essex, CM20 2BU United Kingdom

ABI 3100, 3130XL, 3500, 3500XL, 3730, 3730XL

ABI GeneAmp 9700, Veriti, Proflex, Multiblock system

ABI QuantStudio 5, 6, 7, 12K, ViiA7, 7900 HT FAST

ABI Veriti 384-well

Agilent SureCycler 8800

Analytik Jena Flexcycler2, T1 Thermocycler, Tadvanced, TProfessional, Trobot, Qtower

Bioer Gene Touch 384

Biometra Tadvanced

Bio-Rad C1000, C1000 Touch, S1000

Bio-Rad CFX384, CFX384 Touch

Corbett Research (Qiagen) Palm Cycler 384

Eppendorf MasterCycler ep 384 / Pro 384

GE Healthcare Amersham MegaBACE 4000

Idaho Technology LightScanner

MJ Research (Bio-Rad) DNA Engine, Tetrad, Tetrad 2, Dyad, Dyad Disciple

MJ Research (Bio-Rad) PTC200, 220, 221, 225, 240

peqLab peqSTAR 384X

Sensoquest Labcycler

Techne PCRmax Alpha cycler 1, 2, 4

Techne Prime, PrimeG, Prime Elite, Prime Elite Satellite

Techne TC412, TC512, Genius, Genius Quad, Flexigene

NOTE: We have compiled this list based on plate design and cross-referencing compatibility information available from Eppendorf, Bio-Rad, 4titude, and other sources. Each instrument listed may have variations in block design and configuration. Compatibility is not guaranteed; this is meant as a general guide.

All Instrument brand names, and model names listed above are the trademarks of their respective owners.

For any queries about this guide please contact:

Email: support@3crbio.com

Tel: 01279 940 983

www.3crbio.com